GenBody COVID-19 IgM/IgG

Rapid differential detection kit for IgM and IgG against COVID-19 in human serum, plasma and whole blood

2020.03.13 (Rev.0)

INTENDED USE

GenBody COVID-19 IgM/IgG device is a chromatographic immunoassay kit for the rapid and differential detection of immunoglobulin M (IgM) and immunoglobulin G (IgG) against COVID-19 using serum, plasma and whole blood.

EXPLANATION OF THE TEST

The recombinant COVID-19 antigen was coated on the membrane and anti-human IgM and IgG monoclonal antibody was conjugated the gold particles, respectively. When the specimen existing anti-COVID-19 antibodies is loaded into a sample well (S), the antibodies are complexed with anti-human IgM (or IgG) gold conjugate. And this complex migrates and captured by the immobilized recombinant COVID-19 antigens to make a visible band in the test line regions, M and G. The solution continues to migrate to the control line (C) region that binds a control conjugate, thereby producing another red line. GenBody COVID-19 IgM/IgG can detect the antibodies against COVID-19, so that the device is suitable for the diagnosis of COVID-19 infections.

MATERIALS PROVIDED

GenBody COVID-19 IgM/IgG kit contains the following components:

1. Test device individually foil-pouched with a desiccant
2. Assay solution in dropping bottle
3. Capillary tube for sample loading
4. Instructions for Use

MATERIALS REQUIRED BUT NOT PROVIDED

1. Specimen collection container
2. Micropipette
3. Disposable pipette tips
4. Lancets (for finger prick whole blood only)
5. Centrifuge (for plasma only)
6. Watch or timer

PRECAUTIONS


1. The presence of humidity may decrease the stability of the reagents. Thus, please carry out the test immediately after removing the device from the foil pouch.
2. Do not use the kit after the expiration date. Do not freeze the kit.
3. For in vitro diagnostic use only. Do not re-use the test device.
4. Wear protective gloves while handling samples and wash hands thoroughly after the test.
5. Dispose gloves, swabs, test tubes, and the used strips properly after the test, in accordance with GLP.
6. Do not eat or smoke while handling specimens.
7. Decontaminate and dispose of all specimens, in a biohazard container.

SPECIMEN COLLECTION AND STORAGE

1. Specimen to be tested should be obtained and handled by standard methods for their collections.
2. Serum: Allow the blood to clot, then centrifuge to separate the serum.
3. Plasma: Collect the whole blood into the tube containing anticoagulants such as heparin, citrate, or EDTA. Centrifuge the blood and separate the plasma.
4. Whole blood: whole blood should be collected over heparin, citrate, or EDTA. Mix the blood by inversion and use it to the test. If fingertip blood is used to the test, prick the finger and collect the blood by a capillary tube. And then, load the blood onto the sample well (S) of the test device.
5. All specimens should be tested as soon as early they are prepared. If necessary, they may be stored at 2-8°C for up to 24 hours or at -20°C for longer periods.


TEST PROCEDURE

1. Place all specimens, test devices, and assay solution at room temperature prior to testing (15~30min).
2. [Capillary tube use] Using a capillary tube, add 10 µl of serum/plasma or load 20 µl of whole blood up to black line into the sample well (S).
[Micropipette use] Add 10 µl of serum/plasma or load 20 µl of whole blood into the sample well (S) directly. Add 3 drops (approx. 100 µl) of assay solution into the buffer well (S) in the device.
3. After 10 minutes, interpret the test results.
Please do not read the results after 15 minutes of this testing.

INTERPRETATION OF THE RESULTS

1. Negative: ONLY one band in the control line (C). No COVID-19-specific IgM and IgG were detected. Re-test in 3-5 days if COVID-19 is suspected.
2. IgM Positive: two bands appear in the test line (T) and control line (C) in the left side of device.
3. IgG Positive: two bands appear in the test line (T) and control line (C) in the right side of device.
4. IgG and IgM Positive: each two bands appear in the test line (T) and control line (C) in both side of device.
5. Invalid result: If at 15 minutes, the red band does not appear in the control line (C), the result is considered invalid. If the test is invalid, a new test should be performed with a new patient specimen and a new test device.

[Use quantitative analyzer]


- Using Confiscope G20 is optional.
- Please refer to instructions for Use in the analyzer package.

Confiscope G20

STORAGE & EXPIRATION

1. GenBody COVID-19 IgM/IgG kit should be stored between 2 to 30 °C (35.6 to 86 °F).
2. Expiration date of this kit is 24 months after its manufacture date.

LIMITATIONS OF THE TEST

GenBody COVID-19 IgM/IgG is designed for the detection of IgM and IgG antibodies against novel coronavirus (2019-nCoV). This kit can provide a fast and simple results but, do not completely exclude the possibilities of false positive or false negative results caused by various factors. For confirmation, please make a final decision with clinical symptoms, other testing results, and doctor's assessment, collectively.

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GenBody Inc.
MT Promedt Consulting GmbH Altenhofstr. 80